



Run Info

Host Name	GXB02243 (localhost)
Experiment Name	GZFN211103
Sample ID	UL_DNA
Run ID	0caf19ef-7d95-443e-bff9-324a83ceec97
Flow Cell Id	FAO61000
Start Time	September 29, 10:30
Run Length	3d 0h 3m

Run Summary

Reads Generated	581.6 K
Passed Bases	10.56 Gb
Failed Bases	1.14 Gb
Estimated Bases	12.05 Gb


Run Parameters

Flow Cell Type	FLO-MIN106
Kit	SQK-RAD004
Initial Bias Voltage	-180 mV
FAST5 Output	Enabled
FASTQ Output	Enabled
BAM Output	Disabled
Active Channel Selection	Enabled
Basecalling	on
Specified Run Length	72 hours
FAST5 Reads per File	4000
FAST5 Output Options	zlib_compress,fastq,raw
FASTQ Reads per File	4000
Mux Scan Period	6 hours
Reserved Pores	0 %
Basecall Model	High-accuracy basecalling
Read Filtering	min_qscore=7


Versions

MinKNOW	20.06.9
MinKNOW Core	4.0.3
Bream	6.0.7
Guppy	4.0.11

Cumulative Output Reads




Cumulative Output Bases




Read Length Histogram Estimated Bases - Outliers Discarded

Estimated N50: 38.23 Kb




Read Length Histogram Basecalled Bases - Outliers Discarded

Estimated N50: 37.76 Kb




Read Length Histogram Estimated Bases

Estimated N50: 38.23 Kb




Read Length Histogram Basecalled Bases


Estimated N50: 37.76 Kb




Duty Time Grouped




Duty time Categorised




Mux Scan Grouped




Mux Scan Categorised




Temperature History.




Bias Voltage History.




Translocation Speed



QScore



Disk Write Performance



Run Debug Messages

- The sequencing run has finished, but basecalling may continue October 2, 10:33
- Mux scan for flow cell FAO61000 has found a total of 89 pores. 82 pores available for immediate sequencing October 2, 04:43
- Performing Mux Scan October 2, 04:41
- Mux scan for flow cell FAO61000 has found a total of 162 pores. 133 pores available for immediate sequencing October 1, 22:41
- Performing Mux Scan October 1, 22:38
- Mux scan for flow cell FAO61000 has found a total of 352 pores. 232 pores available for immediate sequencing October 1, 16:37
- Performing Mux Scan October 1, 16:35
- Mux scan for flow cell FAO61000 has found a total of 836 pores. 439 pores available for immediate sequencing October 1, 10:34
- Performing Mux Scan October 1, 10:31
- Mux scan for flow cell FAO61000 has found a total of 130 pores. 106 pores available for immediate sequencing October 1, 05:12
- Performing Mux Scan October 1, 05:09
- Mux scan for flow cell FAO61000 has found a total of 197 pores. 152 pores available for immediate sequencing September 30, 23:09
- Performing Mux Scan September 30, 23:06
- Mux scan for flow cell FAO61000 has found a total of 460 pores. 282 pores available for immediate sequencing September 30, 17:05
- Performing Mux Scan September 30, 17:03
- Mux scan for flow cell FAO61000 has found a total of 1010 pores. 476 pores available for immediate sequencing September 30, 11:02
- Performing Mux Scan September 30, 10:59
- Mux scan for flow cell FAO61000 has found a total of 165 pores. 129 pores available for immediate sequencing September 30, 04:10
- Performing Mux Scan September 30, 04:08
- Mux scan for flow cell FAO61000 has found a total of 346 pores. 233 pores available for immediate sequencing September 29, 22:07
- Performing Mux Scan September 29, 22:04
- Mux scan for flow cell FAO61000 has found a total of 713 pores. 385 pores available for immediate sequencing September 29, 16:03
- Performing Mux Scan September 29, 16:01
- Mux scan for flow cell FAO61000 has found a total of 1436 pores. 499 pores available for immediate sequencing September 29, 10:36
- Performing Mux Scan September 29, 10:33
- Starting sequencing procedure September 29, 10:33
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C September 29, 10:30
- Disk /data has 2479 GB space remaining September 29, 10:30